

# LAPAROSCOPIC CIRCULAR STAPLERS



QUALITY CARE THROUGH INNOVATION



## Laparoscopic Circular Staplers

A FULL RANGE OF ADVANCED CIRCULAR STAPLERS TO COVER ALL SURGEONS NEEDS

Standard  
&  
Extra long



Product Code	Description	Diameter	Staple Height	Closed Staple Height Range	Pack Qty.
FCSLWAF25	Extra Long Circular Stapler with Tilt Top Anvil	25mm	4.8mm	1.6-2.4mm	1/box
FCSLWAF29	Extra Long Circular Stapler with Tilt Top Anvil	29mm	4.8mm	1.6-2.4mm	1/box
FCSLWAF32	Extra Long Circular Stapler with Tilt Top Anvil	32mm	5.0mm	1.8-2.6mm	1/box

Product Code	Description	Diameter	Staple Height	Closed Staple Height Range	Pack Qty.
FCSME21	Circular Stapler with Tilt Top Anvil	21mm	4.5mm	1.4-2.2mm	1/box
FCSME24	Circular Stapler with Tilt Top Anvil	24mm	4.5mm	1.4-2.2mm	1/box
FCSME26	Circular Stapler with Tilt Top Anvil	26mm	4.8mm	1.6-2.4mm	1/box
FCSME29	Circular Stapler with Tilt Top Anvil	29mm	4.8mm	1.6-2.4mm	1/box
FCSME32	Circular Stapler with Tilt Top Anvil	32mm	5.0mm	1.8-2.6mm	1/box
FCSLWBE34	Circular Stapler without Tilt Top Anvil	34mm	5.0mm	1.8-2.6mm	1/box



Panther Healthcare Medical Equipment Co., Ltd  
Beijing | New York | São Paulo | Toronto | Mumbai

E-mail: [orders@pantherhealthcarecanada.com](mailto:orders@pantherhealthcarecanada.com)  
Tel: 1-833-293-8971  
45 Midpark Crescent, London, ON N6N 1A9  
[pantherhealthcarecanada.com](http://pantherhealthcarecanada.com)

[www.pantherhealthcarecanada.com](http://www.pantherhealthcarecanada.com)

# LAPAROSCOPIC CIRCULAR STAPLERS



**Unique 90 degrees Tilt Top Anvil**  
For easy & safe removal through anastomosis,  
minimizing potential tissue trauma




**Improved cutting washer**  
Provides surgeon with audible  
and tactile feedback



**Cartridge and Anvil Self-Lock**  
Prevents tissue from slipping

**Hermetic design**  
Hermetic design  
for laparoscopic surgery,  
prevents the loss  
of pneumoperitoneum



**Dual-Speed Function**

There is less rotation required during the final stages of closure.  
This provides improved instrument control.



**Range of Closure Indicator**

**Dual safety — Patent design**  
Prevents  
accidental misfiring