

HIGH VOLUME HEMORRHOIDAL CIRCULAR STAPLER



QUALITY CARE THROUGH INNOVATION

| Product | Code | Description |
|---------|-----------|---|
| | FCSSWBE31 | Disposable Hemorrhoidal Circular Stapler 31mm, dual safety |
| | FCSSWBE32 | Disposable Hemorrhoidal Circular Stapler 32mm, dual safety |
| | FCSSWBE33 | Disposable Hemorrhoidal Circular Stapler 33mm, dual safety |
| | FCSSWBE34 | Disposable Hemorrhoidal Circular Stapler 34mm, dual safety |
| | FCSSWAE31 | Disposable Hemorrhoidal Circular Stapler 31mm, dual safety, with transparent cartridge |
| | FCSSWAE32 | Disposable Hemorrhoidal Circular Stapler 32mm, dual safety, with transparent cartridge |
| | FCSSWAE33 | Disposable Hemorrhoidal Circular Stapler 33mm, dual safety, with transparent cartridge |
| | FCSSWAE34 | Disposable Hemorrhoidal Circular Stapler 34mm, dual safety, with transparent cartridge |
| | FCSSME32 | High Volume Disposable Hemorrhoidal Circular Stapler 32mm, dual safety, With transparent cartridge and sound feedback |
| | FCSSME34 | High Volume Disposable Hemorrhoidal Circular Stapler 34mm, dual safety, with transparent cartridge and sound feedback |

High Volume HEMORRHOIDAL CIRCULAR STAPLER

Specifically studied for the treatment of hemorrhoidal prolapse and for the S. T. A. R. R. technique




Panther Healthcare Medical Equipment Co., Ltd
Beijing | New York | São Paulo | Toronto | Mumbai

E-mail: orders@pantherhealthcarecanada.com
Tel: 1-833-293-8971
45 Midpark Crescent, London, ON N6N 1A9
pantherhealthcarecanada.com

www.pantherhealthcarecanada.com

HIGH VOLUME HEMORRHOIDAL CIRCULAR STAPLER



Patent design:
Transparent cartridge,
enhanced visibility and
control over captured tissue



2 big casing ports
offer extra visibility
when pulling tissue
into the cartridge



Enlarged indicator window,
convenient for precise
adjustment of staples
closure height



Closing knob: easy
to open and close,
convenient for
operation

Anvil

The biggest case
in the market for
a 34 hemorrhoidal
stapler: 27cc



Markings, indicates
the distance
of anastomotic stoma
from the anal verge

Panther
Healthcare

FCSSME 32

External safety

Firing
handle

Sound feedback, indicates
the stapler is ready
to be fired,
safer operation

Patent design: Dual safety,
avoids misoperation,
much safer



Obturator



Purse-string
Suture Anoscope



Suture
Threader



Circular Anal
Dilator



Butterfly Anal
Dilator