



Open Linear Stapler



Available in four different lengths - 30mm, 45mm, 60mm and 90mm



Two different staple heights: 3.5mm and 4.8mm available, suitable for different tissue thickness



Manual or automatic placement of the alignment/tissue retaining pin



Lockable handle - the handle locks in the back position after firing confirming completed firing cycle



Narrow shaft allows for greater access



Tissue cutting guide for safe and convenient cut line

Colour Indicator Window

Different colours highlight different positions within the firing process



Primary status



Semi-closed status



Fully-closed status



Post-firing status

Product Codes

OPEN LINEAR STAPLER



Product Code	Description	Staple Length	Staple Height
FLSLF30-3.5	Linear Stapler	30mm	3.5mm
FLSLF30-4.8	Linear Stapler	30mm	4.8mm
FLSLF45-3.5	Linear Stapler	45mm	3.5mm
FLSLF45-4.8	Linear Stapler	45mm	4.8mm
FLSLF60-3.5	Linear Stapler	60mm	3.5mm
FLSLF60-4.8	Linear Stapler	60mm	4.8mm
FLSLF90-3.5	Linear Stapler	90mm	3.5mm
FLSLF90-4.8	Linear Stapler	90mm	4.8mm

CARTRIDGES FOR LINEAR STAPLER



Product Code	Staple Length	Staple Height
FLSLH3.5A	30mm	3.5mm
FLSLH4.8A	30mm	4.8mm
FLSLH3.5B	45mm	3.5mm
FLSLH4.8B	45mm	4.8mm
FLSLH3.5C	60mm	3.5mm
FLSLH4.8C	60mm	4.8mm
FLSLH3.5E	90mm	3.5mm
FLSLH4.8E	90mm	4.8mm